2020 Air in circuit

Permission to print: Yes

Incident type Good Catch No Harm Incident

Type of incident: Management Category Air in circuit

Description: HLM using Inspire 6 primed for use. Lines clamped by registrar ready to be used. Arterial

cannulation commenced; and line primed with Plasmalyte by rolling fluid up before and during connection of the cannula to the arterial line. Incidental finding of a air bubble (unknown origin) inside the arterial line (outside of the surgeons vision). Arterial line clamped to patient straight away and line was de-aired using a syringe onto the 22fr EOPA cannula. Line traced back to pump and arterial filter check for further air bubblesnone found. Queried if air bubble introduced while connecting the cannula to the arterial

While there is a checklist item that the circuit is de-aired, a more detailed check of the

line.

GOOD CATCH - what we observers noticed the air bubble in the arterial line and the line was clamped off to the

patient straight away

Preventive actions Thorough check through the lines by all parties pre and post cannulation

better arterial line before cannulation should have been done.

Manufacturer advised: No

What could we do

Discussed with team: Yes

Ext Authority Advised No

Hospital incident filed: No

Knowledge issue No

Rule issue No

Skill issue No